

C20orf166-chromosome 20 open reading frame 166 Gene

Reference: PTXBC110595

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C20orf166; C20orf166; MIR133A2HG; dJ353C17.1; uncharacterized protein MIR1-1HG; MIR1-1 host gene protein; MIR133A2 host; MIR1-1 host gene		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 166
GENE ID	128826
ORIGIN SPECIES	Human
ACCESSIONS	BC110595
VECTOR	pDONR223
NCBI SYMBOL	C20orf166
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagctgtagctgtcattgatggcgcgtgcggccggcgcaggtccgcagcgtggagaggaccagcaggtggcgccgggagagcccggtcggcagctggtcagctccaagtcagccctgaaatgtccattactcacaagaaaaggaaaatgctcatttgaagaa
attcttctctttgtcaatgctgaggcattctctcagccccagccccagtgctccagtggtgaagggaacagctgacggggaagttcagcacatctgtctgaccagggtggagcgatgccagtcctgagctgggaacgcctctctctatgctggagccactgctga