



Reference: PTXBC110597

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LGICZ; LGICZ1; ZAC; ZAC1; zinc-activated ligand-gated ion channel; ligand-gated ion channel subunit; ligand-gated ion channel zinc-activated 1; ligand-gated ion-channel receptor L2; zinc activated ligand-gated ion channel 1; zinc activated ion channel		

For research use only.

GENE DESCRIPTION	zinc activated ligand-gated ion channel
GENE ID	353174
ORIGIN SPECIES	Human
ACCESSIONS	BC110597
VECTOR	pDONR223
NCBI SYMBOL	ZACN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]