

KIAA1602-KIAA1602 Gene

Reference: PTXBC110599

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1602; Cep169; FP1193; nck-associated protein 5-like; NCKAP5-like; centrosomal protein of 169 kDa; NCK associated protein 5 like		

For research use only.

GENE DESCRIPTION	KIAA1602
GENE ID	57701
ORIGIN SPECIES	Human
ACCESSIONS	BC110599
VECTOR	pDONR223
NCBI SYMBOL	KIAA1602
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccagccagctgggggtcctggaacccaaggccaggagagggatgatggcagcatggagccaggcactgcccaggagcttctgcaccgactgcccggagctggaggcagagaactggcacttgcaggccaacgaaaccagcgggagacttatgagcgtgtctg
gacgaggtctgtgggtctgtagtgaggactgggggatgtggtcatctgctcctggcagaagctgggtcagctgatggctctgctcgggctttctgtcccagggttccaacatgtgtacaggcgtgtgctgaaccagaagggtgagtaa