

EPN3-epsin 3 Gene

Reference: PTXBC110605

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	epsin-3; EPS-15-interacting protein 3; epsin 3		
For research use only.			

GENE DESCRIPTION	epsin 3
GENE ID	55040
ORIGIN SPECIES	Human
ACCESSIONS	BC110605
VECTOR	pDONR223
NCBI SYMBOL	EPN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacacctcgcactccggcggcagggtgaagaacatcgtgcacaactactcaggcagaaatcaaggtgcgcgagggccaccagcaatgacccctggggccccctagtctcgtcatgtccgagatcgtgacctgaccttcaacacagtggccttcaccgaagtcagggcatgctgtggcggcggtcaatgacagcggcaagaactggcggcacgtgtacaaggcttaacattgctggactactctcaagacgggtctccgagcgggtggccccaccagtccgagagaacctctacaccatccagacactcaaggacttccagtacatcgaacggcagggcaaggaccaggcgtcaacgtgcgcgagaaggtcaagcagggtgatggccctgtcaaggatgaggagcgggtgcgcgagggcgaacccacgcctcaagaccaaggagcgcattggcactggaggcatcggcattggcagttgggcagctgggcttcagcccgctacggcgaggactacagccgtctccgggggtccccgtctcctctacaactctctctctctcgttcaccccgctatacctcgcacgtggagcagggccggcctcagacgtcaggggaaggaactgcagctgcagctggccctgcctcagcctgagccgtgaggagcagagaagcctgtccccccagcctccacaggacgagagacctgcagctgcagctgctctgcgcctgagccggcaggagcagagagaagggtgagggtctctggcagggtgatggctccccatggccaatggtgcaggggccgtgtgtccaccatcagcgggacagagagcctgagagagaagagaagagaggagggagaagc taaaaaccagccagggccccccactctgagccgggtgcctctccccccacagtcctccatctggacttggctga