



Reference: PTXBC110615

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD209 molecule; CD209 antigen; CDSIGN; CLEC4L; DC-SIGN; DC-SIGN1; C-type lectin domain family 4 member L; HIV gp120-binding protein; dendritic cell-specific ICAM-3-grabbing non-integrin 1; dendritic cell-specific intercellular adhesion molecule-3-grabbing		

GENE DESCRIPTION	CD209 molecule
GENE ID	30835
ORIGIN SPECIES	Human
ACCESSIONS	BC110615
VECTOR	pDONR223
NCBI SYMBOL	CD209
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]