

KRTAP9-2-keratin associated protein 9-2 Gene

Reference: PTXBC110619

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP9.2; KRTAP9.2; keratin-associated protein 9-2; keratin-associated protein 9.2; ultrahigh sulfur keratin-associated protein 9.2; keratin associated protein 9-2		

For research use only.

GENE DESCRIPTION	keratin associated protein 9-2
GENE ID	83899
ORIGIN SPECIES	Human
ACCESSIONS	BC110619
VECTOR	pDONR223
NCBI SYMBOL	KRTAP9-2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccactgtgtcccttctgtcagcctacatgtctcaggaccacctgtcaggaccacctgtggaagcccaccactgtgaccacctgcagcagcacatcctgtgccagcccgctgtgtgtccagctgtccagcctgtgccgccaactcctgtcaaaacacgtg
ctgtaggaccacctgtgccagccacctgtgtgaccagctgtgccagcctcctgtcagcagcacacctgtgccagccacctgtgtgggtccagctgtgtgccaaccagctgtgggtccagctgtggccagagcagctcctgtgcacctgtgtactgcagaagaacctgtac
taccacgactgtgtcctgctgtgtgtgcctaaaccagagctgtgtgctccaactgtgccagccgtgtgtgccagccacctgtgtgagaccacctgtgcaggaccactgtgtccagccacctgtgtgtccagctgtgccagccttctgtgtgtgta