



Reference: PTXBC110621

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	heparan-sulfate 6-O-sulfotransferase 2; HS6ST-2; heparan sulfate 6-O-sulfotransferase 2		

For research use only.

GENE DESCRIPTION	heparan sulfate 6-O-sulfotransferase 2
GENE ID	90161
ORIGIN SPECIES	Human
ACCESSIONS	BC110621
VECTOR	pDONR223
NCBI SYMBOL	HS6ST2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]