

PAGE5-P antigen family, member 5 (prostate associated) Gene

Reference: PTXBC110803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT16; CT16.1; CT16.2; GAGEE1; PAGE-5; P antigen family member 5; P antigen family, member 5 (prostate associated); cancer/testis antigen 16.1; cancer/testis antigen family 16, member 1; cancer/testis antigen family 16, member 2; g antigen family E member		

For research use only.

GENE DESCRIPTION	P antigen family, member 5 (prostate associated)
GENE ID	90737
ORIGIN SPECIES	Human
ACCESSIONS	BC110803
VECTOR	pDONR223
NCBI SYMBOL	PAGE5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgcacatgaacaagatcccaatcctcagaagaggaaatgaccaagagctctccagccagttggacgtgtgattgtccagcagcccactgaggaaaaacgtcaagaagagggaaccaccaactgataatcagggtattgcacctagtggggagatcaaaaatgaaggagcacc
tgctgttcaaggagctgatgtggaagctttcaacagggaactggctctgcttaagatagaggatgcacctggagatggtcctgatgtcaggagggggactctgccacttttgatcccaactaaagtgtggaagcaggtgaagggaactatag