

TSSK4-testis-specific serine kinase 4 Gene

Reference: PTXBC111088

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf20; STK22E; TSK-4; TSK4; testis-specific serine/threonine-protein kinase 4; serine/threonine kinase 22E; serine/threonine-protein kinase 22E; testis specific serine kinase 4		

For research use only.

GENE DESCRIPTION	testis-specific serine kinase 4
GENE ID	283629
ORIGIN SPECIES	Human
ACCESSIONS	BC111088
VECTOR	pDONR223
NCBI SYMBOL	TSSK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagtcttgcggcacaagctacatcaactctatcggccattgagagcacatctcgagtatacatcattctggaactggctcagggtggtgatctcctgaatggatccagcgctacgggctgctctgagcccttgcctggcaagtgttctccagctgacctgggcattgcctacctgcacagcaagagcatcgtgcaccgggacttaagtggagaacctgttgcctggacaagtgggagaatgtgaagatatcagacttggccttggccaagatggtgccttctaaccagcctgtgggtgttagcccttcttacgcccaagtgaactgctttccacctcagccagacttactgtggcagcttgccttacgcttcccagagatctacgagcctggccctacaaccttctctgtctgacacctggagcatggcgctcatcctttacactctagtggctgcccatctgccccttgatgacaccaatctcaaaaagctgctaagagagactcagaaggaggtcactttccagctaacataccatctccaggagtgcaagaacctgatcctcagatgctacgccaagccactaagcgtgccaccattctggacatcatcaaggattcctgggtgctcaagttccagcctgagcaaccccatgagatcaggctgcttgaggccatgtgccagctccacaacaccactaaacagcaccaatccttgcaattacgacctga