

**39509-membrane-associated ring finger (C3HC4) 2 Gene**

Reference: PTXBC111388

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

**For research use only.**

---

<b>GENE DESCRIPTION</b>	membrane-associated ring finger (C3HC4) 2
<b>GENE ID</b>	51257
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC111388
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	39509
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgacgacgggtgactgctccacccccggctccctgtgtgactgctcggcagccctgctctccaaggtcgtggaggtacggcctcggacgccccagtagtggcacaggtgactcaagggatggccgctcctccaccgtcctcctgacacaccgagtgat  
ggctctttctgcccgatctgcatgagggagcgaacgggagtgcttctgctccccctgtggctgcaccggcagcctgggtgcccgtgcataagagctgtctggagaagtgcttctcctcctaacaccagctactgcgagctgtgccacacggagttgcatggagaacggcctcga  
ccccacagagtgctgaagaccggggccgacgagaagcggacactgtctgacatggtgtttcctgctcaccaccgctggcccatctcaggctggttggcctgcgcccggaccaccctcggctccacagccagctggagccgtgggtctcattg  
ccctcaccatgccccttccacctatgtctctggacgtgctcctcctccctaccactgccagctgtactccgagtgagagaagaccaaccagaaagtgcctgaagatccggaggcggacagccccaggccccccagcattctccactggcagctggactctgaagaagg  
tgcagaggagacaccagtatga