

MUCL1-mucin-like 1 Gene

Reference: PTXBC111421

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SBEM; mucin-like protein 1; mucin like 1		

For research use only.

GENE DESCRIPTION	mucin-like 1
GENE ID	118430
ORIGIN SPECIES	Human
ACCESSIONS	BC111421
VECTOR	pDONR223
NCBI SYMBOL	MUCL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttcttagcagtcctgtactctggagtttccatcttctgtctctgcccagaatccgacaacagctgctccagctgacacgtatccagctactggctctgctgatgatgaagccctgatgctgaaaccactgctgcaaccactgcgaccactgctgctctaccactgcaacc
accgctgcttaccactgctctaaagacattccagttttacccaaatgggttgggatctcccgaatggtagagtggtccctga