

BARX2-BARX homeobox 2 Gene

Reference: PTXBC111432

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein BarH-like 2; BarH-like homeobox 2; BARX homeobox 2		
For research use only.			

GENE DESCRIPTION	BARX homeobox 2
GENE ID	8538
ORIGIN SPECIES	Human
ACCESSIONS	BC111432
VECTOR	pDONR223
NCBI SYMBOL	BARX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcactgccacgcccagctgagctgagctgcccggcagctcaagcagccagcggcgctacaagactttcatgacgacgagatcctccaaggagacctgcgattactttgagaaacttccctctactcgtgtgccgtgctgtgctgctgcgaccaagccccctgattctgtacgggctccccctcctgcccgcataccgctcctcctcgggtgatcaccggccagcccactgtcatctcccactggtccctgccaccccgggaatgccaggcactgtctgccaccaggtcaccgaggcgggtctctgctgagggccccagggggcaggccctagccagcagcga gtcagagacggaacgcccagccccacagaagaagccccgggagtcgacccatcttcaccgagctgcagctcatgggctggagaagaattccagaagcagaagtattgtcaaccccagacaggttgacttggtcagctctctggactcactcagctgcaggtgaagacc tgggtatcagaatcgaggatgaaatggaagaaatggttcttaaggtggacaggaagcaccacaaaacccaaggtcggcccaagaagaactcctcccatcagaagagattgaagctgaagagaagatgaacagccaggccaggggtcaggagcagctggagccctctca ggggcaggaggagctctgtgaagcacaggaaacgaacgtgatgtcccttagagatggcagagccaccagacccggccaggagttgccaataccctcttgcgaacccccaccattaaagctaa