

C21orf123-chromosome 21 open reading frame 123 Gene

Reference: PTXBC111539

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf123; NCRNA00175; PRED80; COL18A1 antisense RNA 1		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 123
GENE ID	378832
ORIGIN SPECIES	Human
ACCESSIONS	BC111539
VECTOR	pDONR223
NCBI SYMBOL	C21orf123
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctttgtgccctgccgggtcctggcacggcgaggtgcgtcaggccagccagcagcaaggcccggtcctggcttctgccatctattctgtcctgcaactctatctggagcctccatccctcaaaagctgcaggaagaggcctgcaggcaagagaaggagctgca
ggagtgcgaatttaagaatggaaggccgaggttttccaatctaattaacagcatgacaagccaggatgggtctgacccgctggagccttgacagccagaaacgcggtgtcctgtccgctcctgagacgcagggctgtcactgcaacacagccccaacttccacacggg
cggactctggctcctcaggcaccgcagtgccctctcctccctccagctctcctga