

LOC284837-hypothetical protein LOC284837 Gene

Reference: PTXBC111558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC284837		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC284837
GENE ID	284837
ORIGIN SPECIES	Human
ACCESSIONS	BC111558
VECTOR	pDONR223
NCBI SYMBOL	LOC284837
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatgcagactctacatcaagggaacaagtcggaggccgaggggtccacggcactccgacccacggggccagggtggcacaacgccccaacgctgcaatcctgggaaggtcacacgcgcctccgggacccacgtgacctcaaggagtgaggagacaca
tcctcgggggtgacgcccctgcagatggagttcccatctgtgccaccttgcaggcctgggaaagctcggggccaggcggcaccgggcaagctgagtgccagtggtgtgtgcaacacacgggcaaatctgaaacctgggacggacgagctgcacccctgggaaaga
gaacagaggaagccgggttatcaacgtcctcctgtgctcagttccagaacaagccagactcagtgctcaggagacccctgtgtgctccggagagcacgggacccgccctttgtcaggtgaggaactgaacccagagcagtgaggactggctgggacccctgggctc
agccctgggccccttccctctgccatctcagacgggcacagccagcggcctcaccacaacggcctccagactggccctcaggggcccccaggcacaatcctctcactttctcctcaagccctgccttcccccttgcactgtgacccaggggctgctagcagctgggctgggc
tccaaagcttgaagctgcttccaggtcttccaggggcgccggatggctgctcactcaggacgccaacacacaggaaggacccatgcacggatctgacttccctcacaggctctgccttcaactaatgtccctgtggggcagctggggctcccggccctccctgcc
acagccgtgaggttgcctctag