



Reference: PTXBC111562

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ECSM4; MRB; roundabout homolog 4; roundabout homolog 4, magic roundabout; roundabout, axon guidance receptor, homolog 4; roundabout guidance receptor 4		

For research use only.

GENE DESCRIPTION	roundabout homolog 4, magic roundabout (Drosophila)
GENE ID	54538
ORIGIN SPECIES	Human
ACCESSIONS	BC111562
VECTOR	pDONR223
NCBI SYMBOL	ROBO4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]