



Reference: PTXBC111564

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf1; C17orf1B; FHOD4; FMNL; KW-13; formin-like protein 1; CLL-associated antigen KW-13; leukocyte formin; formin like 1		

For research use only.

GENE DESCRIPTION	formin-like 1
GENE ID	752
ORIGIN SPECIES	Human
ACCESSIONS	BC111564
VECTOR	pDONR223
NCBI SYMBOL	FMNL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr