



Reference: PTXBC111698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RUTBC2; small G protein signaling modulator 1; RUN and TBC1 domain containing 2; RUN and TBC1 domain-containing protein 2; small G protein signaling modulator 1 protein		

For research use only.

GENE DESCRIPTION	small G protein signaling modulator 1
GENE ID	129049
ORIGIN SPECIES	Human
ACCESSIONS	BC111698
VECTOR	pDONR223
NCBI SYMBOL	SGSM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]