

KIAA0040-KIAA0040 Gene

Reference: PTXBC111699

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA0040		
For research use only.			

GENE DESCRIPTION	KIAA0040
GENE ID	9674
ORIGIN SPECIES	Human
ACCESSIONS	BC111699
VECTOR	pDONR223
NCBI SYMBOL	KIAA0040
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcactatgtccatgtccacagagtaactactcaaccaaggacaacacacagactaagtggtccagtgaggggcagtcaccagcgtggacaattcttgatactgtcttggcagctatgtccaatagcaatgctcctactgcagaccaggcatgcctccacctgtctctgg
cataccacacatgcaagcacaagaacattatccatacatcctaattggtcccaagtggtgacatgcacgtacacacacacacacaaattcaggtacaggtacgtgggcaagtattctgctcatcaaatggcattggctatgtactttgtcaggggaagtacattatctacagtca
caaaaatgtctcatgggaaagccttgcagattcagacacatatatacaatttctaaccagcaaggcccccatacaccatctattccataa