



Antibodies and Proteins

**ProteoGenix SAS**  
Espace Européen de l'Entreprise  
15 rue de La Haye  
67300 Schiltigheim, France  
(0)3 90 20 54 70 - Fax +33 (0)3 88  
contact@proteogenix.fr  
www.ProteoGenix.science

### SCAND3-SCAN domain containing 3 Gene

Reference: PTXBC111742

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SCAND3; Buster4; ZFP38-L; ZNF305P2; ZNF452; dJ1186N24.3; SCAN domain-containing protein 3; SCAN domain containing 3; transposon-derived Buster4 transposase-like protein; zinc finger BED domain-containing protein 9; zinc finger protein 305 pseudogene 2; zi |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | SCAN domain containing 3                                                                          |
| <b>GENE ID</b>          | 114821                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC111742                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | SCAND3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)

Page 1/1