



Reference: PTXBC111765

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf47; CCDC35; DHCD1; DNHD1L; dynein heavy chain domain-containing protein 1; DNHD1 variant protein; coiled-coil domain containing 35; dynein heavy chain domain 1-like protein; dynein heavy chain domain 1		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 47
GENE ID	283294
ORIGIN SPECIES	Human
ACCESSIONS	BC111765
VECTOR	pDONR223
NCBI SYMBOL	C11orf47
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr