

FLJ40448-hypothetical protein FLJ40448 Gene

Reference: PTXBC111769

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ40448
GENE ID	339059
ORIGIN SPECIES	Human
ACCESSIONS	BC111769
VECTOR	pDONR223
NCBI SYMBOL	FLJ40448
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatggcctatgcatgtaccagtcctccatggccaccctggcatttcaggcacctgtgtgcgcaggtcaggccccagaccggcctctctctctgtgcacccgcacctggcaggaagcggcagcaccatctggtgcctggcagcctgtctctctccagctgcctgtggc
cccactggactcagccccagcagggaggccttcggctgcagggcagggcctggcgggcgccacctgtgtgtggggcactcaggagccctggcttcacctgtctctcggcgccacatgtggacagctcaccctccagcaggggtccggaccaagctctgacccaccgggtt
ctgagtcctgacagggcgtgtggacagcttccttccacccggcactgtctggcaggagcttcggggccctcacctggggcagcaactcctgcaccaccagccacggccggccacgagaaagtgcagaacagctcaggccctctcacaactcctggagctcagggacaggtag