

DUS4L-dihydrouridine synthase 4-like (*S. cerevisiae*) Gene

Reference: PTXBC111774

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DUS4; PP35; tRNA-dihydrouridine(20a/20b) synthase [NAD(P)+]-like; protein similar to E.coli yhdg and R. capsulatus nifR3; dihydrouridine synthase 4 like		

For research use only.

GENE DESCRIPTION	dihydrouridine synthase 4-like (<i>S. cerevisiae</i>)
GENE ID	11062
ORIGIN SPECIES	Human
ACCESSIONS	BC111774
VECTOR	pDONR223
NCBI SYMBOL	DUS4L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaaggttatgggcttcttaataaacaagccagagctgttcaagacatgttgaacaagtaagaatcaagtggaaacccctggatttccagtttctataaaataaggatccatgatgaccttaaaagaactgtagatctttgtcaaaaggctgaagcaacaggagtttcatggatt
acagtccatggaagaactgtgaaagaagacatcagccagtgactatgattccattaaaataaaggaaaatagtctatactgttaattgctaatggagacatcagaagctaaaggagcagaaaatgtgtggcggattactgggacagatgggtgtgatgttgcagaggactcttag
caaaccggccatgtttgctgatatgaggaaacccactgaaatgcactctgggactgggttgacattgctctgaactcgggactcctacatgtgttccatcaacatttaattacatgatggaaaagataactcaaggcaggaaaaagggtatttaattgctctgtcaagcacatcagcaa
tcatagattacctacagaccattatggcatttga