

UNC13B-unc-13 homolog B (C. elegans) Gene

Reference: PTXBC111781

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MUNC13; UNC13; Unc13h2; protein unc-13 homolog B; homolog of rat Munc13 (diacylglycerol-binding); munc13-2; unc-13-like; unc-13 homolog B		

For research use only.

GENE DESCRIPTION	unc-13 homolog B (C. elegans)
GENE ID	10497
ORIGIN SPECIES	Human
ACCESSIONS	BC111781
VECTOR	pDONR223
NCBI SYMBOL	UNC13B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr