

LOC348262-hypothetical protein LOC348262 Gene

Reference: PTXBC111795

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC348262
GENE ID	348262
ORIGIN SPECIES	Human
ACCESSIONS	BC111795
VECTOR	pDONR223
NCBI SYMBOL	LOC348262
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatgggagaagagccacgggagagggcaggctggaggcagggtaccctaaagcagcagtcggcagtggtgagagccagcaggggcaggcaggggcctggccagctctacctgttacctgagctctgcttctctgtaacgggagctccagcaggcagcatgtccct
gtgggacctcaacccaacaggcctttccctccgtgctctggtgttctgtggcctggagtcctggcaggaggactggcagagagacccagagtcgaagaaggagaggtgactttgtgagcaaaactgggtgctgcccgtgggtgggagccctggcccttttggacctcaactct
ggcctgggaggggacagagttccaggcgtgggagctggtttctgctcttggctgtttggcccttgagccgtgggagttcttatcacgtggtgtttgaggcctggaccattgacatgagcgggaatgagccagagaggactgaagcctcagtgctcctggccctctgtga