

## KIF1A-kinesin family member 1A Gene

Reference: PTXBC111799

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | HCA2; HEL-S-7a; SAC; SACI; Sacy; hsAC; adenylate cyclase type 10; 3',5'-cyclic AMP synthetase; AH-related protein; ATP pyrophosphate-lyase; epididymis secretory sperm binding protein Li 7a; germ cell soluble adenylyl cyclase; testicular soluble adenylyl c |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | kinesin family member 1A                                                                          |
| <b>GENE ID</b>          | 547                                                                                               |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC111799                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | KIF1A                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]