

TAC4-tachykinin 4 (hemokinin) Gene

Reference: PTXBC111858

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HK-1; HK1; PPT-C; tachykinin-4; endokinin; preprotachykinin-C; tachykinin 4 (hemokinin)		

For research use only.

GENE DESCRIPTION	tachykinin 4 (hemokinin)
GENE ID	255061
ORIGIN SPECIES	Human
ACCESSIONS	BC111858
VECTOR	pDONR223
NCBI SYMBOL	TAC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccttgcctgcctgccttctctgatggagctgtccgtgtgcactgtggcaggtgatggaggaggaacagacactcagcactgaagcagagacctgggaaggcgtggccccagcattcagctccagctgcaggaggtgaagacgggcaaggcaagccagttcttgggctgatggggaagcgaagtgggaggtgagtgga