

MYH7-myosin, heavy chain 7, cardiac muscle, beta Gene

Reference: PTXBC112171

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMD1S; CMH1; MPD1; MYHCB; SPMD; SPMM; myosin-7; cardiac muscle myosin heavy chain 7 beta; myHC-beta; myhc-slow; myopathy, distal 1; myosin 7; myosin heavy chain beta-subunit; myosin heavy chain slow isoform; myosin heavy chain, cardiac muscle beta isoform		

For research use only.

GENE DESCRIPTION	myosin, heavy chain 7, cardiac muscle, beta
GENE ID	4625
ORIGIN SPECIES	Human
ACCESSIONS	BC112171
VECTOR	pDONR223
NCBI SYMBOL	MYH7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr