



Reference: PTXBC112203

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR167; MGRF; MRGE; mas-related G-protein coupled receptor member E; G-protein coupled receptor 167; MAS-related GPR, member E; Mas-related G protein-coupled receptor F; mas-related G protein-coupled MRGE; MAS related GPR family member E		

For research use only.

GENE DESCRIPTION	MAS-related GPR, member E
GENE ID	116534
ORIGIN SPECIES	Human
ACCESSIONS	BC112203
VECTOR	pDONR223
NCBI SYMBOL	MRGPRE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]