

FLJ14816-hypothetical protein FLJ14816 Gene

Reference: PTXBC112205

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC84931		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ14816
GENE ID	84931
ORIGIN SPECIES	Human
ACCESSIONS	BC112205
VECTOR	pDONR223
NCBI SYMBOL	FLJ14816
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacggaacggaggccccggctggggcaggcgcctgctgatacaggaaggacaaaggctcgacgattccggcttctcttagctgctctctgcaagctcttctccccatgggtctctgatgagccaacaagacgggcagggagcccgatctcacagatggg
gaagctgaggccagcgagccttccccagcacagaagagagctcctggcattccagcaggggggactggaagaaaaccaggacctggcttaaccacctcactcttcagatggggaatagctggagatacaggaatggcagagtagattctgactggaggttccccagggtg
ccccctgctggttgggggaacacaggagccacgacctccccaggttggctgcatgtttatggattccagaatatttaacagctgagtggaatggcactggctcaaggccgcctgttccctcccatctaccacaaggtatcttctgggaagcagggaagacagtgtag