

C6orf54-chromosome 6 open reading frame 54 Gene

Reference: PTXBC112215

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf54; HGC6.1.1; NCRNA00300; KIF25 antisense RNA 1		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 54
GENE ID	26236
ORIGIN SPECIES	Human
ACCESSIONS	BC112215
VECTOR	pDONR223
NCBI SYMBOL	C6orf54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagacatttttaaaataatggagtcctccagggcggtgcgtgctgtggcctgtgtcctcaccgctttggaccaggtgcgtgtctccaccagaccaaggactcacagagctgcctgggggacctggccccacagtcacccgtccgtcatcagcctcaatgccatcag
cccgatgaatgctgtccatcgtcgtcagcagccaaggaggtctctcctctaaagcacctggagcccttggatggfagcaggtgatgtggcgatgagtggtgcctcgtggggccttggacgagcgtggacaccccaatccccagactggccattgccgtggagcagcgtgag
tgtgacatggagctcggtgtcgtgctgcagaggacgctcgcgctgtacgctcatgatagccaggacccctccacatgtcacttagccaaggatgcagcccgccctggggctttctacccagcccaggacctgcagggacaaggacccctcagcgagatgcagctctcat
gaggcctga