

LIMS3-LIM and senescent cell antigen-like domains 3 Gene

Reference: PTXBC112233

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PINCH-3; LIM and senescent cell antigen-like-containing domain protein 3; LIM and senescent cell antigen-like domains 3; LIM-type zinc finger domains 3; particularly interesting new Cys-His protein 3; pinch 2; LIM zinc finger domain containing 3		

For research use only.

GENE DESCRIPTION	LIM and senescent cell antigen-like domains 3
GENE ID	96626
ORIGIN SPECIES	Human
ACCESSIONS	BC112233
VECTOR	pDONR223
NCBI SYMBOL	LIMS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccttctcaggccgagcgcgccctgcattatccagagaacgaagaatccccgagcagcccttaacactgtccacgaggccaatgggaccgaggacgagagggtgtttcctaaactgcagcgcaggcagtgacgtgaaagtctacaaggagttctgtgacttttatgcgaa
attcaacatggccaacgccctggccagcgcacttgcgagcgtgcaaggcggtttgcgccgctgagacgaicgtgaacagtaattgggagctgtaccatgagcagtggtttcgtgtgcgtcagtgcttcagcagtgccagaaggactttctatgaggaacgaacgtga