

C13orf30-chromosome 13 open reading frame 30 Gene

Reference: PTXBC112237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf30; protein FAM216B; family with sequence similarity 216 member B		
For research use only.			

GENE DESCRIPTION	chromosome 13 open reading frame 30
GENE ID	144809
ORIGIN SPECIES	Human
ACCESSIONS	BC112237
VECTOR	pDONR223
NCBI SYMBOL	C13orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacaaaactggaaaagacaacaaaagctttggaatgtccacaactcttttattcaggttcctccatctatgacactcttactaaaggccctcaaccaagggaacagcgtactttacagcattatgaggatttacaactccaggccccagtgaggagccctgcagaccg
ctacattcacagccttcagcaccaacagctgcttggtatattactcaacgggaagccttgcttctgtacttagagattcaaccaagagagcctcagccaaggtagctcctcaagaacattccccgaaaacttcagccatgacaagaagatgccatcagtaactacgtgtatctgt
ggttctacctaggcccaagttaaaggcgccaagtgtcaggaactga