

ATOX1-ATX1 antioxidant protein 1 homolog (yeast) Gene

Reference: PTXBC112248

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	copper transport protein ATOX1; ATX1; HAH1; ATX1 antioxidant protein 1 homolog; metal transport protein ATX1; antioxidant 1 copper chaperone		

For research use only.

GENE DESCRIPTION	ATX1 antioxidant protein 1 homolog (yeast)
GENE ID	475
ORIGIN SPECIES	Human
ACCESSIONS	BC112248
VECTOR	pDONR223
NCBI SYMBOL	ATOX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgaagcagcaggttctctgtggacatgacgtgtggaggctgtgtgaagctgtctctcgggtcctcaataagctggaggagttaatgacattgacctgcccaacaagaaggtctgcatgaaatctgagcacagcatggacactctgttgcaaccctgaagaaaacaggaaagac
tgtttcctaccttgccctgagtag