

CCL7-chemokine (C-C motif) ligand 7 Gene

Reference: PTXBC112258

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FIC; MARC; MCP-3; MCP3; NC28; SCYA6; SCYA7; C-C motif chemokine 7; chemokine (C-C motif) ligand 7; monocyte chemoattractant protein 3; monocyte chemotactic protein 3; small inducible cytokine A7 (monocyte chemotactic protein 3); small-inducible cytokine A		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 7
GENE ID	6354
ORIGIN SPECIES	Human
ACCESSIONS	BC112258
VECTOR	pDONR223
NCBI SYMBOL	CCL7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagcctctgcagcacttctgtgtctgtctcacagcagctgtttcagccccagggcgttgcagccagttgggattaatacttcaactacgtctgtctacagattatcaataagaaaatccctaagcagaggctggagagctacagaaggaccaccagtagccactgtcccggaagctgtaattcaagaccaaactggacaaggagatctgtgtgacccacacagaagtgggtccaggacttatagaagcacctggacaagaaaacccaactcgaagcttga