

GCHFR-GTP cyclohydrolase I feedback regulator Gene

Reference: PTXBC112262

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GFRP; HsT16933; P35; GTP cyclohydrolase 1 feedback regulatory protein; GTP cyclohydrolase I feedback regulator		

For research use only.

GENE DESCRIPTION	GTP cyclohydrolase I feedback regulator
GENE ID	2644
ORIGIN SPECIES	Human
ACCESSIONS	BC112262
VECTOR	pDONR223
NCBI SYMBOL	GCHFR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctacctgctcagcaccagatccgatggaggtgggccccactatgtggggcgaacagtcggatccagagctgatgcagcatctggggcttcaagagaagagccttgggaaacaactttatgaatactacgtcgtgacccctcccccatagctcctggacaagctg
gaacgcaggggcttcctgtgtctgagcatgacgggggtgggccagacgctgggtgtgtgtctgcacaaggagtg