



Reference: PTXBC112265

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B4Gal-T5; BETA4-GALT-IV; beta4Gal-T5; beta4GalT-V; gt-V; beta-1,4-galactosyltransferase 5; UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 5; UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 5; UDP-galactose:beta-N-acetylglucosamine beta-		

For research use only.

GENE DESCRIPTION	UDP-Gal:betaGlcNAc beta 1,4- galactosyltransferase, polypeptide 5
GENE ID	9334
ORIGIN SPECIES	Human
ACCESSIONS	BC112265
VECTOR	pDONR223
NCBI SYMBOL	B4GALT5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]