



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

FAM19A3-family with sequence similarity 19 (chemokine (C-C motif)-like), member A3
Gene

Reference: PTXBC112268

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM19A3; TAFA-3; TAFA3; chemokine-like protein TAFA-3; family with sequence similarity 19 (chemokine (C-C motif)-like), member A3; family with sequence similarity 19 member A3, C-C motif chemokine like		

For research use only.

GENE DESCRIPTION	family with sequence similarity 19 (chemokine (C-C motif)-like), member A3
GENE ID	284467
ORIGIN SPECIES	Human
ACCESSIONS	BC112268
VECTOR	pDONR223
NCBI SYMBOL	FAM19A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgagagggctcgagcggaaactggagcacggcgctgcgtctggcactgtgcctggcctgctgtggaccacctggaccttggctgccttcgagcctccactctgcacatgtctgtgtgcagcaggcgacctgcgagggtgattgcggctcacctgtctgcaaccggaaaccgc
atcgagagggcctccagacgttgaaatgctctcttcttggcagctgtggccgacacggcggaagccctctgcctgtggcagactctgttgcctcctctgtctgtgagacacactagatgctcactccgctctctctcccaacacgctctcccaacgctctcctgcagagatggtg
ctctgatagatgaacccctcctccctccggagagatgttaagctgtctcccgacgtctcgatgagatctcagacagctgtggcagacaaatgcaaaacacaaatgcacatagatctcttgggtgctcgcgcctctgcagaaagaagctctgactga

More info: Please contact us for any question: contact@proteogenix.fr