

PIN4-protein (peptidylprolyl cis/trans isomerase) NIMA-interacting, 4 (parvulin) Gene

Reference: PTXBC112281

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rotamase PIN4; peptidyl-prolyl cis-trans isomerase Pin4; PPIase PIN4; EPVH; PAR14; PAR17; peptidyl-prolyl cis-trans isomerase NIMA-interacting 4; eukaryotic parvulin homolog; hEPVH; hPar14; hPar17; parvulin; parvulin-14; parvulin-17; peptidyl-prolyl cis/t		

For research use only.

GENE DESCRIPTION	protein (peptidylprolyl cis/trans isomerase) NIMA-interacting, 4 (parvulin)
GENE ID	5303
ORIGIN SPECIES	Human
ACCESSIONS	BC112281
VECTOR	pDONR223
NCBI SYMBOL	PIN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatggcgggcttctaaagggctgtacggcaactggagcagttcaggttcaacaacaagcttcaagatgccgccaaggaaaaagtggtcttgaaaaagcggggaagggggagcagcctctgggagtgacagtgctgacaagaaggctcaaggtcccaaggtggt
ggcaatgcagtaagggtcagacacattctatgtgaaaaacatggcaaaatcatggagccatggaaaagttaaagctctgggatgagattcaatgaagtgccgcacagtatatgagataaagccaggcaagggggtgacttgggttgatgaccagagggtccatggtgggaccatt
caagaagcagcatgttccttgcctgttaagtggttgataagcctgtgtttacagaccaccggttaagacaaaattggatatcatattattatggtcgaagggaagaaataa