

SLC35E2-solute carrier family 35, member E2 Gene

Reference: PTXBC112289

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	solute carrier family 35 member E2		
For research use only.			

GENE DESCRIPTION	solute carrier family 35, member E2
GENE ID	9906
ORIGIN SPECIES	Human
ACCESSIONS	BC112289
VECTOR	pDONR223
NCBI SYMBOL	SLC35E2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtcctcgtgaaacccagcactggaagagctggtcctcgcgaagagaagccgaaaggcaggtcgcctcagctgggctcctgtttggtcaccgaagtgagaagattgtttgccaaagcgcagggcgccacagatgagaacgtactgaccgtcaccatcacgga
gaccacgggtcagtcagacttgggtgtgtggagctcgcggcgctgctacctcagctgtggttcttctcagcttctgcagctccttctcaacaagtacatcctgtcctcgtgggagcgagcccgatgtaggtgcgggtgcagatgctgtccaccacggtatcgggtgtgtg
aaaaccctcgttcttctgtttatcagcacaaggccggccttctacccaccaacttcttatgacgatgctgtttggtgctgatgaggtttgcaactgtgtttggtttggtcagcctgaaaaatgtggcggttctgtgtgagacggtgaagagctccgcccatcttcacggtg
atcatgtcggatgattctggggagttacacagagctccagtgatcgggagggcgggaagagctcagctacaaccaggacgtggtgctgctccttgacagacgagtcctgttccactcagagcgtcagggcgtacgcccctatgggaaaaatctcccggtgacttcag
gtccgcagggcctgcaccagtcgcttgcctatggcgtcgcctgccccatgcttctcgtcgtcga