

IFNA1-interferon, alpha 1 Gene

Reference: PTXBC112300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFL; IFN; IFN-ALPHA; IFN-alphaD; IFNA13; IFNA@; interferon alpha-1/13; IFN-alpha 1b; IFN-alpha-1/13; interferon alpha 1b; interferon alpha-D; interferon-alpha1; leIF D; interferon alpha 1		

For research use only.

GENE DESCRIPTION	interferon, alpha 1
GENE ID	3439
ORIGIN SPECIES	Human
ACCESSIONS	BC112300
VECTOR	pDONR223
NCBI SYMBOL	IFNA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctcgccctttgcttactgatggctctggtggtcagctgcaagtcaagctgctctggtggtgatctccctgagaccacagcctggataacaggaggaccttgatgctctggcacaatgagcagaatctctctctctgctgatggacagacatgactttgattcccc
aggaggagtttgatggcaaccagttccagaaggctccagccatctctgctccatgagctgatccagcagatctcaacctttaccacaaaagattcatctgctctggtggtgaggacctctagacaaattctgcaccgaactctaccagcagctgaatgacttggaagcctgtgtgatg
caggaggagagggtgggagaaactccctgatgaatgcggactccatcttggtctgtgaagaatacttccgaagaatcactctctatctgacagagaagaatacagcccttgctcctgggaggtgtcagagcagaatcatgagatccctctttatcaacaacttgaagaagatta
aggaggaaggaataa