

MLN-motilin Gene

Reference: PTXBC112314

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	promotilin; prepromotilin		
For research use only.			

GENE DESCRIPTION	motilin
GENE ID	4295
ORIGIN SPECIES	Human
ACCESSIONS	BC112314
VECTOR	pDONR223
NCBI SYMBOL	MLN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtatcccgtaaggctgtggctgctctgctggtgcatgcagctgcatgctggcctccagacggaagccttcgtcccatcttcacctatggcgaactccagaggatgcaggaaaaggaaacggaataaaggcgaagaaatccctgagtgatggcagaggtctggggagga
aggctctgtagaccctgcggagcccatcagggagaagaagaacgaatgatcaagctgactgctcctctggaattggaatgaggatgaactccagacagctggaaaagtaccggccaccctggaaggctgctgagtgagatgcttcccagcatgcagccaagtga