



Reference: PTXBC112336

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DUTT1; SAX3; roundabout homolog 1; deleted in U twenty twenty; roundabout, axon guidance receptor, homolog 1; roundabout guidance receptor 1		

For research use only.

GENE DESCRIPTION	roundabout, axon guidance receptor, homolog 1 (Drosophila)
GENE ID	6091
ORIGIN SPECIES	Human
ACCESSIONS	BC112336
VECTOR	pDONR223
NCBI SYMBOL	ROBO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1