

SPIN4-spindlin family, member 4 Gene

Reference: PTXBC112360

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TDRD28; spindlin-4; spindlin family member 4		
For research use only.			

GENE DESCRIPTION	spindlin family, member 4
GENE ID	139886
ORIGIN SPECIES	Human
ACCESSIONS	BC112360
VECTOR	pDONR223
NCBI SYMBOL	SPIN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctccaaccgtgctccgatggggtagatggcgtgtccgcatacctgatgaagaaaggcacacccacagggaagcaacggcgcaagcccactttctcactcgttagaacatcgtgggctgccgcatcacaacggcgtggaagggaagcaacgagccagtgagcagtggaagggtactgtctcgagcaggtttccgtgaagcccactctttacatcattaaatgatggcaagatagtggtatggactagaactgcacccgataagagagtttagcgtctagagatcctctgagagagtgccaactcctcgtatcgattcacgactggcagattccctgattggcaaggcagtgagcatgtgttgaagggaacatggtaccaaggatgaatggaagggtatggtcctggcgcgagctcctgtgatggatactgtgtttacatcacctacgagaagatcctgttctctatatgtacacgctgcttgatgactacaagatgggtgacttacgcattatccagattccaaactatttccctacagcagaacaggagcctggagaggtggtcgacagtctctgtgggcaagcaggtggagcatgccaaagatgacgggtccaagagaactggcattttatacatcaagtggcgcaagccatctgtttactcattaaagtttgatgatattcacattatgtctatggttggtgaaaactccttaa