

SERP1-stress-associated endoplasmic reticulum protein 1 Gene

Reference: PTXBC112364

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	stress-associated endoplasmic reticulum protein 1; ribosome associated membrane protein 4; ribosome-attached membrane protein 4; stress associated endoplasmic reticulum protein 1		

For research use only.

GENE DESCRIPTION	stress-associated endoplasmic reticulum protein 1
GENE ID	27230
ORIGIN SPECIES	Human
ACCESSIONS	BC112364
VECTOR	pDONR223
NCBI SYMBOL	SERP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcgccaagcaaggatcgatggccaacgagaagcagcaagaacatcaccagcgcggaacgtcgccaagacctgagaaatgccccgaagagaaggcgtctgtaggacctggtattggctctctcattttgtgtctgtggtctgcaattttccagattattcaagta
tcaggatggcatgtga