

C6orf12-chromosome 6 open reading frame 12 Gene

Reference: PTXBC112395

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf12; NCRNA00171; ZNRD1-AS; ZNRD1AS; ZNRD1 antisense RNA 1		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 12
GENE ID	80862
ORIGIN SPECIES	Human
ACCESSIONS	BC112395
VECTOR	pDONR223
NCBI SYMBOL	C6orf12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagaaaacccagagcgtagaaggactaaacagctactgtttgtacctgagagagatgaaacacatagaaggcacatacatcaaacaggaaaggctggagaatttaagaacaaaccatttagacaagtactacaaccccccaatgaaacaaaattacaaaaattatg
ccagaaggcgtgcatccagaacgcacagagaagaagcaggtaaatgagagggaacagatgcagactaaggatcaccagggaacgatgattcgaggagagaactgcagaacaaagactcaaggagagaatcttgagaagaagccagagccagctactacatagagaagc
atgagagagtaaaagagataaaggagtttgaagagttatgctgtatcttttccaacctgcagtagaagtcggattaaagtgaagtattcttggataagtcagaatggagagaagggtgaatacaattgtgaaaccatacaagaaaattcttgcaatgccacccttttaagaagtca
aataggaaaaataagagattaa