

KGFLP1-keratinocyte growth factor-like protein 1 Gene

Reference: PTXBC112977

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	keratinocyte growth factor-like protein 1
GENE ID	387628
ORIGIN SPECIES	Human
ACCESSIONS	BC112977
VECTOR	pDONR223
NCBI SYMBOL	KGFLP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgattcgggtccagattctggcaggagggagtttgggaltcagatctgaaaaagcactagactggaagaggacgcgatggagtcggagccgtggcgggtacaaaaaccagagccggggaagcgcgggtggggaggaagcagcggatggactttacctgcgcac
gcgtcagccatctccgcgacagtggtggccaccgcgactggtcgtgaagtgtcggcgtgtccggcgtccgctgggacccgggtgctcgcctgagctcagcttctctatctgtacgggtggacaagtacagtaaccctgccctcaagacgggccagggctgtgagcag
gggtccacgcttagagcaggcacctattgtgcaggccctgagatgggtctgactcagttcctgcggggaacttcaccagtgaccagtcagtgccctcagttaa