

FLJ37201-hypothetical protein FLJ37201 Gene

Reference: PTXBC112983

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ37201
GENE ID	283011
ORIGIN SPECIES	Human
ACCESSIONS	BC112983
VECTOR	pDONR223
NCBI SYMBOL	FLJ37201
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaaatagaaaacattgcatagagccaaaatgaagacctgaccctgtgttgattgaatgggttcgacaacaagaagtaaaatagacaccactgactggtttgttgctgtaacaagctagaatataccatgaagaactgaactttgaaggagagtgtgaatattcaggaggctggctgcagaaatttaaaaagtgtcttggaatcaagtatcttcaaatctgtgtgaaaaggcttctctgatgaaagctactgaaaattacattgatgaatttgcgaaggaataatctgatgaaatcttagccctgaacaaatttacagtgctgatgaacagctttgtactggtgctatgttcctagaaaactttaacaatggctaaacaaagagcaacaacagatgtcaaaagataactaagcaaaagttaactgttctgtgagtgctaatgctgcaggcacacacaagataaaattggcagtgattggaaatagcctataccaagggtgtttaaagggtgtgcataatttacctgtacattgttatgcacacaagaaagcatcagtaacgggaatctttccaactggttga