

LOC148413-hypothetical LOC148413 Gene

Reference: PTXBC112998

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC148413		

For research use only.

GENE DESCRIPTION	hypothetical LOC148413
GENE ID	148413
ORIGIN SPECIES	Human
ACCESSIONS	BC112998
VECTOR	pDONR223
NCBI SYMBOL	LOC148413
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtccctagcggagcgccgggtgccgtgctgagagcgcctgctgtgcgccccgagcggggctgggactctccaagatgccacgttcgacagagacccggatcgccgaagctcgctctcgaaaggcctcggtctcacgccctgcccgctcctgggttcacgggttttca
tcacctgcggctgctctcgatcgaccacagctgtgcaggaggcagggatctgtgtcgtcagttaccggaaacctttgccagtagtacaggaccacggcgtgtagctcaggatgtctcgtctgtgagttacagctgcacgctctccaggaaagaaggaaatttctctctctg
gaaacccaccacacagctggtttctcattggctgcttgccttccctga