



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF151-ring finger protein 151 Gene

Reference: PTXBC113014

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

RING finger protein 151

For research use only.

GENE DESCRIPTION	ring finger protein 151
GENE ID	146310
ORIGIN SPECIES	Human
ACCESSIONS	BC113014
VECTOR	pDONR223
NCBI SYMBOL	RNF151
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctggcggggtatgacatcaacctctctgccagccctcctgacagcaactctgtgtctgctccgtgccatgggggttcaagaggccgaaggttgccatcgagccacatctctgcataaaagtcatctccgggtgctagccagacaaaagacctgtcgtgctgtaggaaagaggtg
gaaaagggaagaaaggtgtgcacatgaataaacctccgaaacaacttggccctctggaagtcgaagtcgaagaacccgacgtgctgctcatlgtgacatgcccttggcccatcgcaaggggacaccaggactatgcccttltgactaacggcctgcccacaggagggtgtcacctc
cgaggtgtccgtctgggaccttggcagagaccgcgcagcattgccgaagaaggctccagcagcgtgtccccctgggtctcggggccacctgtgaccgcccgacgctgctgccacaactgttaccgggagctgtcaccaagcgttgagcgtgtccgagagcgcctgtcgccgc
ctgtgctgtccctctgtggcgtgtgtgctgtgctggaccacaccacatgtcgttcgtagagagctgtggcgagctcgaacactctggaagaaacaccgctctgtggagggtgtccccacagaggaggccgagctgtcccccagagggaacgttggggctgaggtgtgtggg
ggagcccgaggtccaactcctgtgaataag

More info: Please contact us for any question: contact@proteogenix.fr