

LOC285733-hypothetical LOC285733 Gene

Reference: PTXBC113027

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC285733
GENE ID	285733
ORIGIN SPECIES	Human
ACCESSIONS	BC113027
VECTOR	pDONR223
NCBI SYMBOL	LOC285733
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgagcaaaagccggagccccagaagaaacaagtacagactcagaggaaagctgccctggtcctgagtggtgactccatggtcccggtgggtctgtgtgtgttgcaatgagctctgtgctgctcagcttcatgaggagctccctggctggttctgttcttgggtcttctcccc
atgactttgctgctcctcctcatgcctacttcaggatcaaaactgattgaggttaatgaagaactgtcccagaactgtgatgcacaataatcccaaggatggctcttccctgtaccagagaatgaaatggactga